



PIUMA

NANOINDENTER



Optics11 B.V.
De Boelelaan 1081
1081 HV Amsterdam
The Netherlands

For more information,
please visit our website:
www.optics11.com

or contact us at:
info@optics11.com

Office telephone:
+31 (0)20 598 79 17

Piuma Nanoindenter Application note:

Measuring the mechanical properties of cartilage immersed in PBS

The Piuma Nanoindenter is specifically developed to non-destructively characterize the mechanical properties of soft materials, such as tissues, biomaterials, hydrogels and cells. Using a novel proprietary optical fiber-top technology, the Piuma provides great force sensitivity, accuracy and precision while remaining very easy to use.

Application note AN-1503, written by Ernst BreeI, Optics11 BV
Published on 3 November 2015 on www.optics11.com

Cartilage is a widely studied material, as defects due to trauma or disease can severely impact quality of life. At the basis of developing new therapies lies a thorough understanding of both biological and mechanical properties of cartilage. This note describes the mechanical characterization of rat cartilage using the Piuma Nanoindenter.

Introduction

To characterize the mechanical properties of rat cartilage, two experiments are performed: varying the probe displacement from 7 to 20 micrometer, resulting in an increasing depth of indentation, and mapping of a section of the cartilage to explore local variations in stiffness.


Methods and materials

The Piuma Nanoindenter is fitted with a 7 N/m probe with a 88 μm radius extended spherical tip for the first experiment and a 7.95N/m probe with a 24.5 μm radius extended sperical tip for the second experiment. The displacement profile is varied from 7 to 20 μm from the surface, with 2 $\mu\text{m}/\text{s}$ loading and un-loading displacement rates and 6 s holding time.

Applying the Hertzian model describing the indentation of elastic bodies, the data in the loading section of the load-displacement curve is used to determine the Youngs Modulus, using a fit of all datapoints from the contact point to 60% of the maximum load point [1].

For these experiments, the Piuma Nanoindenter is placed on top of a regular lab bench; no special stabilization or dampening is used during the experiments.

Sample preparation




Immersed cartilage sample mounted on Piuma X-Y stage.

Bone samples from adult rat are harvested, cut below the metaphyseal compartment, parallel to the cartilage surface, and frozen at -80°C for one week. The sample is thawed at room temperature

and the epiphyseal section is attached to a glass substrate using a thin cyanoacrylate film (see picture above). The sample is fully immersed in PBS (pH 7.4) during the measurement.


Results

To investigate the E-modulus response to varying indentation depth, one series of indentations using 7 to 20 μm displacement profiles are set. The displacement rate and hold time are kept constant.



As the data shows, the resulting indentation depths vary between 1.2 to 5.0 μm . The E-modulus increases with increasing indentation depth and load.

Additionally the spatial E-distribution of a site on the cartilage is examined, using a 12 μm displacement profile.



The scan reveals that the cartilage surface has regions with softer and stiffer properties, when using a 20 μm pitch in X and Y direction.

[1] Hertz, H. R., 1882, Ueber die Beruehrung elastischer Koeper , in Gesammelte Werke (Collected Works), Vol. 1, Leipzig, Germany, (1895)