


Overview

The MicroSquisher has been designed to perform tension and compression testing at low forces. Properties that can be measured when a **spherical aggregate of cells is compressed between parallel plates** (Fig. 1) include **cell viscosity** and **surface and interfacial tensions**. These properties are important because they determine how cells and tissues self-organize during embryogenesis, wound healing, regenerative medicine, and cancer metastases. The compression test provides a controlled mechanical environment in which the aggregate and its cells are deformed. By monitoring the geometry of the aggregate (Fig. 1B) and the time course (Fig. 2) of the force the aggregate exerts on the compression plates, one can determine the effective viscosity of the cytoplasm in the cells, the strength of the interfacial



tension that acts along each cell-cell boundary (Fig. 3) and the strength of the surface tension present along the cell-medium interface³.

The Compression Test

In a typical compression test, a spherical aggregate of cells is placed between parallel plates. The plates are then moved together under and the aggregate is compressed. In this case, the compression is typically carried out relatively rapidly (over several seconds) and, as a result, individual cells in the aggregate deform as the aggregate deforms.





Fig. 1 – The compression process. Immediately following compression, the aggregate and its cells are visibly deformed. After annealing, the cells become nearly isotropic in shape again, but the aggregate remains compressed (from Brodland et al., 2009).

Cell Aggregate (Spheroid) Compression Testing

With time, the cells reshape and reorganize within the aggregate and the platen force decreases. [Specific parts of the force-time curve can be associated with particular cell properties^{2,3,5}](#). Specifically, the effective viscosity of the cells (μ) can be determined from F_μ , the rapid decrease in force from B1 to B2. The force decay F_{cc} from B2 to C is caused by reshaping of the cells inside the aggregate, and it can be used to calculate the cell-cell interfacial tension γ_{cc} . The steady-state force F_{cm} provides data for calculating the cell-medium and cell-platen interfacial tensions γ_{cm} and γ_{cp} , respectively. These surface and interfacial tensions arise from the action of various sub-cellular mechanical components including acto-myosin systems, membrane contractions and cell-cell adhesion systems^{1,4}, and the measured forces can be used to quantify their mechanical activity.



1. Brodland GW (2002) DOI: 10.1115/1.1449491
2. Brodland GW. (2003) *Biorehology* 40, 273-277
3. Brodland GW et al. (2009) DOI: 10.2976/1.3175812
4. Lecuit and Lenne (2007) DOI: 10.1038/nrm2222
5. Yang and Brodland (2009) DOI 10.1007/s10439-009-9649-6

©Copyright 2016, CellScale Biomaterials Testing



[CellScale Biomaterials Testing](#) is the industry leader for precision biomaterial and mechanobiology test systems. Our products are being used at world-class academic and commercial organizations in over 30 countries around the globe.

Our [mechanical test systems](#) allow researchers to characterize the mechanical properties of biomaterials. Our [mechanobiology technologies](#) provide insights into the response of cells to mechanical stimulation.

[CellScale's technologies](#) are improving human health by helping researchers discover the causes of disease, improve medical treatments and devices, and advance regenerative medicine and other basic science research.

[Visit our website](#) or [contact us](#) to learn how our innovative products can help you achieve your research and development goals.